

Biopsy forceps

The entire product range of HSO



Product properties

HSO disposable biopsy forceps are the **inexpensive** and **safe alternative** to reusable biopsy forceps. The **different jaw parts** allow a **gentle and efficient extraction** of meaningful tissue samples.

HSO disposable biopsy forceps are manufactured according to European standard norms. The use of **Metal Injection Moulding Technology** – a metal powder injection moulding process – guarantees the **easy movement** and **precision** of the movable jaw parts.

Thanks to the **high quality workmanship**, HSO disposable biopsy forceps can be introduced easily and gently for tissue sampling.

Coated and uncoated metal spirals are at your disposal for the sample collection in the specialist areas > of bronchoscopy, colonoscopy and gastroscopy.

- Reliable
- Ease of movement
- Precise



ORDERING OPTIONS GASTROSCOPY

Head	Coating	Diameter	Length	Shipping unit	Order No.
With mandrel, not serrated	None	2.3 mm	1.8 m	100 pieces	HS02318S
Without mandrel, not serrated	None	2.3 mm	1.8 m	100 pieces	HS02318
With mandrel, not serrated	PE-coated	2.3 mm	1.8 m	100 pieces	HS02318SC
Without mandrel, not serrated	PE-coated	2.3 mm	1.8 m	100 pieces	HS02318C
Without mandrel, serrated	None	2.3 mm	1.8 m	100 pieces	HS02318T
Without mandrel, serrated	PE-coated	2.3 mm	1.8 m	100 pieces	HS02318CT
With mandrel, serrated	PE-coated	2.3 mm	1.8 m	100 pieces	HS02318SCT

ORDERING OPTIONS COLONOSCOPY

Head	Coating	Diameter	Length	Shipping unit	Order No.
With mandrel, not serrated	None	2.3 mm	2.3 m	100 pieces	HS02323S
Without mandrel, not serrated	None	2.3 mm	2.3 m	100 pieces	HS02323
With mandrel, not serrated	PE-coated	2.3 mm	2.3 m	100 pieces	HS02323SC
Without mandrel, not serrated	PE-coated	2.3 mm	2.3 m	100 pieces	HS02323C
Without mandrel, serrated	None	2.3 mm	2.3 m	100 pieces	HS02323T
Without mandrel, serrated	PE-coated	2.3 mm	2.3 m	100 pieces	HS02323CT
With mandrel, serrated	PE-coated	2.3 mm	2.3 m	100 pieces	HS02323SCT

ORDERING OPTIONS BRONCHOSCOPY

Head	Coating	Diameter	Length	Shipping unit	Order No.
Without mandrel, not serrated	None	1.8 mm	1.8 m	00 pieces	HS01818
Without mandrel, not serrated	None	1.8 mm	1.2 m	100 pieces	HS01812
Without mandrel, not serrated	PE-coated	1.8 mm	1.8 m	100 pieces	HS01818C
Without mandrel, not serrated	PE-coated	1.8 mm	1.2 m	100 pieces	HS01812C

CE 0408