

# Cleaning brushes

## Disposable cleaning brushes



### Product characteristics

**HSO cleaning brushes** have been developed for the efficient cleaning of test tubes, cannulas, nozzles, endoscopes and other medical devices.

The benefits of **HSO cleaning brushes** at a glance:


- Single use guarantees maximum cleaning effect
- Gentle bristle tips prevent damage to working channels etc.
- A flexible pulling tube and unique positioning of the bristles allow simple, efficient forward and backward movements
- A secure grip and adhesion of the brushes is guaranteed by the welding to the pulling tube – no bonding
- Welded sheathings prevent the entry of fluids into the pulling tube
- Easy handling
- Latex-free

The HSO cleaning brushes are disposable products.

The single use of the HSO cleaning brushes ensures a maximum cleaning effect!

Please always use the proper brush diameter for a comfortable application.

The HSO cleaning brushes are available in different sizes and various versions:



Upon request, we also offer extremely thin brushes for small-caliber working channels and other versions in different sizes for different working channels and working channel openings.

Contact us!

ORDERING OPTIONS HSO CLEANING BRUSHES					
Product	Brush diameter	Brush length	Overall length	Shipping unit	Order No.
1-sided	5.5 mm	20 mm	2,300 mm	100 pieces	HSO18230
1-sided	10 mm	20 mm	110 mm	100 pieces	HSO10040
2-sided	3 mm	20 mm	2,000 mm	100 pieces	HSO17200D
2-sided	5 mm	20 mm	2,300 mm	100 pieces	HSO18230D

