



**OPFix**

**HISO**<sup>®</sup>  
HOSPITAL SUPPLY ORGANISATION



HAND AND FOOT  
RESTRAINT

**SECURE  
IN ALL  
SITUA-  
TIONS**

## Two-piece hand and foot holder

Patients with uncontrolled movements or situations in which such patients must be repositioned often present an enormous challenge for nursing staff and care workers. In these situations, catheters, hoses, etc. often present a great risk of injury to the patient.

By restraining the arms and legs, **OPFix prevents the uncontrolled manipulation** of catheters, drains or life-sustaining equipment for the **patient's own safety**. Thanks to its easy handling, the two-part product set **assists** in the repositioning of patients.

### OPFix HAND AND FOOT RESTRAINT



OPFix **consists of two padded hand and foot holders as well as one restraining strap**. The attachment rings are available in both metal and plastic. For patients with a limited urge to move (e.g. patients with dementia or geriatric patients, in radiology settings or during MRI scans), we recommend the use of the cuffs with the plastic ring. The cuffs with the metal ring should be used on patients with a strong urge to move as well as on patients with a tendency for aggressive behaviour.

#### Order options

Product	Shipping unit	Order No.
OPFix with metal ring	25 pairs	HSOOPFIX1000
OPFix with plastic ring	25 pairs	HSOOPFIX1000P
FixBand	100 pieces	HSOFIXBAND6000

### OPFixSoft HAND AND FOOT RESTRAINT



OPFixSoft noticeably improves comfort for both the patient as well as the practitioner. The **extremely soft material of the cuff protects** the skin against abrasions from sharp edges at the seams and so considerably increases the wearing comfort for the patient.


The **strap with the metal ring is firmly sewn to the cuff**, which makes it easier to apply. It can be tied or knotted depending on the situation and intended use. The product is ready for use with no need for preparatory measures. The cuff is exclusively available with the metal ring.

#### Order options

Product	Shipping unit	Order No.
OPFixSoft	25 pieces	HSOOPFIX800

[www.hso-solutions.at](http://www.hso-solutions.at)

If you have any questions, please contact  
your personal product specialist:



**SALES**

Rudolf Heintel Gesellschaft m.b.H.  
Erdbergstraße 166  
1030 Vienna, Austria

T +43 (1) 403 89 56-0

F +43 (1) 403 89 56-50

[medizintechnik@heintel.at](mailto:medizintechnik@heintel.at)

[www.heintel.at](http://www.heintel.at)



**MANUFACTURER**

MSP Medizintechnik GmbH  
Würtzlerstraße 16  
1030 Vienna, Austria

T +43 (1) 403 89 56-0

F +43 (1) 403 89 56-50

[info@hso-solutions.at](mailto:info@hso-solutions.at)

[www.hso-solutions.at](http://www.hso-solutions.at)

Photo:

© IStock Photo ozgurdonmaz

© IStock Photo Anton Novik

Subject to errors and omissions  
Version 2.0 | Last updated: Nov 2017